

Original Article

Does retensioning of adjustable-loop cortical suspension devices improve performance: A systematic review and meta-analysis

Seth Theodore Campbell¹, Joseph Liu², Anirudh K. Gowd³, Nirav Amin⁴, Matthew Sardelli¹ , Martin Morrison⁵

¹Department of Orthopedic Surgery, McLaren Flint, Flint, Michigan, ²Department of Orthopedics, University of Southern California, Los Angeles, California, ³Department of Orthopedic Surgery, Wake Forest Baptist Medical Center, Winston-Salem, ⁴Premier Orthopedics, Pomona, California, ⁵Department of Orthopedics, Renown Pediatric Orthopedics and Scoliosis, Reno, United States.

ABSTRACT

Objectives: To date, there is conflicting evidence when comparing fixed-loop cortical suspension devices (FLDs) to adjustable-loop devices (ALDs). Some studies indicate that ALDs are inferior to FLD in regard to displacement and failure load while others show that they are biomechanically similar. The purpose of this study is to use a meta-analysis of biomechanical data to compare FLDs to ALDs with and without retensioning. It is hypothesized that retensioning the ALD will allow these devices to be biomechanically equivalent to the FLD in total unloaded displacement and failure load.

Materials and Methods: This study sought to identify all biomechanical studies that compared fixed loops to ALD. A meta-analysis was performed to find the standardized mean difference with retensioning as a covariate.

Results: The analysis of isolated tests showed that retensioning reduced the cyclic ALD displacement in comparison to non-retensioned ALD; however, both the ALD with and without retensioning had significantly higher cyclical displacement and significantly lower failure load compared to the FLD. In the meta-analysis of the animal model data, there was no significant difference between the ALD with retensioning and the FLD.

Conclusion: This analysis suggests that retensioning reduces displacement of an ALD, with displacements measured in animal bone testing showing no significant difference as compared to FLD. However, in the device-only model, the ALD with retensioning and the ALD without retensioning had significantly higher cyclical displacement and significantly lower load to failure compared to the FLD. These data suggest that retensioning may be beneficial. However, there is significant heterogeneity in the pooled studies which limit the strength of this conclusion.

Keywords: Anterior cruciate ligament, Cortical suspension device, Adjustable loop, Biomechanical

INTRODUCTION

In the United States, alone it is estimated that there are approximately 130,000 anterior cruciate ligament (ACL) reconstruction surgeries annually.^[1] There are many techniques used for graft fixation: cortical suspension devices, interference screws, and cross-pin devices, with each having advantages and disadvantages.^[2,3] Due to concerns over the biomechanical strength of the interference screws^[4,5] and the operative complications of the cross-pin devices,^[2,6,7] other methods of fixation, such as cortical suspension devices, have been developed to optimize the biomechanical and insertional qualities.

Two main categories of cortical suspension devices are the fixed-loop device (FLD) and the adjustable-loop device (ALD). FLD can provide biomechanically sound fixation with minimal elongation as the graft incorporates into

the bone. However, these devices also require precise calculations regarding the graft and tunnel length, and they do not permit change in length when encountering differing bone tunnel dimensions.^[8] ALD allows for simpler insertion by avoiding intraoperative calculations as well as minimizing over-drilling,^[9] which maximizes the bone-graft interface. There are concerns, however, that these devices may exhibit increased displacement with cyclical loading when compared to FLD-type devices in some biomechanical studies.^[2,3,9-15]

To date, the literature has provided conflicting evidence when comparing FLDs and ALDs. Some biomechanical studies concluded that the FLD showed less total displacement as well as higher ultimate failure loads.^[2,3,9-14] Other biomechanical studies have concluded that fixed and ALD are biomechanically similar in total displacement as well as failure load.^[8,16-18] One of the unique benefits of an ALD is

*Corresponding author: Seth Theodore Campbell, Department of Orthopedic Surgery, McLaren Flint, Flint, Michigan, United States.
sethcampbell2012@gmail.com

Received: 16 August 2022 Accepted: 16 August 2023 Epub Ahead of Print: 17 April 2024 Published: XXXXX DOI: 10.25259/JASSM_22_2022

This is an open-access article distributed under the terms of the Creative Commons Attribution-Non Commercial-Share Alike 4.0 License, which allows others to remix, transform, and build upon the work non-commercially, as long as the author is credited and the new creations are licensed under the identical terms. ©2024 Published by Scientific Scholar on behalf of Journal of Arthroscopic Surgery and Sports Medicine

the ability to re-tension after intraoperative cyclical loading, which allows the surgeon to retighten the device and remove some of the initial laxity. Few studies, including a recent meta-analysis,^[15] have evaluated the effects of retensioning in their analysis. There is some evidence to suggest that retensioning the ALD could possibly make the ALD biomechanically similar to the FLD,^[8,13,17,18] while a single study concluded that it made no difference.^[13]

The purpose of this study is to compare FLDs to ALDs with and without retensioning to determine whether retensioning the device can recapture some displacement and make it biomechanically equivalent to the FLD. It is hypothesized that retensioning the ALD would allow these devices to be biomechanically equivalent to the FLD in total unloaded displacement and failure load.

MATERIALS AND METHODS

Search criteria and data extraction

PubMed, Embase, and Cochrane Library databases were searched for eligible studies up to July 15, 2020, using the search terms (cortical or suspensory) and “anterior cruciate ligament” and “biomechanical.” Resulting abstracts were reviewed to identify any biomechanical model studies that compared adjustable-loop systems to fixed-looped systems and tested for displacement and or load failure were included for analysis. Any studies that published retractions or studies that were challenged in published letters to the editor were excluded from the study.

Included studies were examined to identify the primary outcomes of interest which were the total displacement (in mm) and failure load (in N) of each device. In cases where the data of interest were provided graphically (rather than in the text or a table), the authors were contacted to request the study data. In these cases where the authors did not respond (two cases) an image analysis tool (ImageJ, NIH) was used to digitally extract the data.^[14,19] In cases where data were reported as median/range, the method of Hoza *et al.*^[20] was used to convert the data into estimated mean and standard deviation. Rylander *et al.* did not evaluate for displacement using methods similar to the other studies and so the displacement data from that study was not included in the analysis.^[21] Study methods were also extracted and documented in tables.

Quality/bias assessment

The methodological quality/risk of bias of the selected articles was assessed independently by two reviewers using a standard checklist and scoring method [Table SI 1].^[22] Briefly, each article was evaluated for five items with each scored as 0 if well addressed, 1 if partially addressed/unclear, and 2 if not addressed. The assessment items focused on (1) sampling methods for the devices included, (2) clarity and consistency in testing protocol, (3) sample size justification, (4) clarity

and consistency in data analysis protocol, and (5) complete reporting of data. A total quality score was calculated and articles scoring $0 < 3$ were assigned a grade low risk of bias, 4–7 rated moderate risk, and 8 or higher rated high risk. Where there was a discrepancy between the reviewer’s ratings, the score was discussed to arrive at a consensus.

Statistical analysis

Data from device-only models (DOMs) were analyzed separately from data gathered using animal models to limit heterogeneity in the testing methods of each study, as demonstrated in Kamelger *et al.*'s^[14] study which showed significant differences between results from these testing modes. Statistical analysis was performed using the open-source software OpenMeta (<http://www.cebm.brown.edu/openmeta/>) and SigmaStat (Systat Software Inc. San Jose California USA).

The primary analysis was a sub-group meta-analysis to examine the within-study differences between ALD and FLD. In this analysis, the subgrouping was based on whether the ADL was retensioned or not. The standardized mean difference and 95th confidence intervals (CIs) for each group were found using a random effects model.^[8] In studies where there were multiple control devices tested, ALDs were matched to FLD using a random assignment process. A secondary meta-analysis grouped all devices into three groups (grouping data across studies): the FLD (Group 1), the ALD without retensioning (Group 2), and the ALD with retensioning (Group 3).

For all meta-analyses, a Wald test was performed to determine if the standardized mean difference for each group (for the within-study difference analysis) or the group means were non-zero. Forest plots were created for each analysis and the I^2 index was used to measure the heterogeneity of included studies.^[23] Index scores $< 25\%$ are considered to indicate low heterogeneity, 20–50% moderate, and $> 75\%$ considered high.^[24] Secondary analyses to identify significant differences between retensioned and not retensioned ADL pooled means from the meta-analysis of within study differences utilized *t*-tests. Analysis of variance analyses with Holm-Sidak *post hoc* testing were utilized in comparing data pooled across studies for Groups 1–3. All analyses were done with a 95% CI.

RESULTS

Using the above search criteria, a total of 274 studies were identified in the search. Among these 33 full-text articles were assessed for eligibility [Figure 1]. Within these 18 studies were excluded after review, leaving 15 final studies included in the study.

The included studies involved a range of devices with the most common control being the ENDOBUTTON™ (Smith and Nephew, Memphis, TN) and the most common adjustable devices being the TightRope® RT (TRT, Arthrex, Naples, FL)




Figure 1: Search strategy using the preferred reporting items for systematic reviews and meta-analyses guidelines. n=Number of studies.

and ToggleLoc™ Device with ZipLoop® Technology (TLZ; Zimmer Biomet, Warsaw, IN) [Table 1].

None of the included studies were determined to present a high risk of bias. Three studies^[9,14,17] were assessed to have moderate potential for bias based on a missing or poorly documented sample size calculation and lack of documentation/potentially non-random sample selection. Please refer to the on line only Supplementary Information (SI), [Table SI 1] for the score for each paper.

Experimental setup

The setup for the “device-only model” studies involved testing the devices in isolation. These studies typically passed the device through a hole in a steel plate rather than bone. After feeding the button through the opening, the button was secured against the plate surface simulating a lateral femoral cortex. The free end of the loop was then put over a rod or a hook. The loop lengths remained equal for both the ALD and FLD in each study but differed between studies.

In the tests utilizing animal bone models (ABMs), bovine tendons were cut to a certain length, then doubled or quadrupled over and sutured together to replicate an ACL construct. Graft measurements were taken to size the bone tunnel. The femoral tunnels were drilled to various distances with a reamer leaving a varying amount of bone bridge at the cortex. The tunnel length is the distance of larger diameter tunnel for the graft. The graft with the attached loop device was pulled through the femoral tunnel, pulled through the smaller diameter tunnel until the button was on the cortex, and the button was then flipped. The ALDs were, then, adjusted so that the graft filled the tunnel. Götschi *et al.*^[25] attached a tendon graft to the device without using animal bone. Smith *et al.*^[8] used both a femur and a tibia in their animal model with the cortical button fixation on the femoral side and an interference screw on the tibia side. However, in their TightRope® RT group, an all-inside technique was used with an ALD on both the femur and tibial side.^[8] Most studies (DOM and ABM) precycled its construct. During this phase, lower forces for lower amounts of

cycles were used to remove any slack from the setup and simulate intraoperative graft cycling. Any low load residual displacement after this was set to zero as the baseline. The total displacement value was recorded after subsequent higher load the cyclical loading. Loading levels and numbers of cycles in each study are shown in [Table 1]. After the cyclical loading protocol, failure testing would begin by pulling the setup by a certain distance per unit time until the setup would fail, this would be recorded as the failure load force.

Five studies retensioned their ALD but only three studies compared a retensioned device to a non-retensioned device. Johnson *et al.*,^[13] Noonan *et al.*,^[17] Singh *et al.*,^[26] and Nye *et al.*^[18] retensioned by pulling on the sutures after precycling. Smith *et al.*^[8] used 200 N of force, or maximal force if 200 N could not be reached, to re-tension the femoral side. Two studies tied a knot on the outside of the adjustable loop.^[10,17] The groups that had a knot but were not retensioned were analyzed in the non-retensioned group.

Meta-analysis of total displacement

The data extracted for the displacement analysis can be found on-line in the SI Data Source Tables ([Table SI 2a] for DOM displacements and [Table SI 2b] for ABM displacements).

In the DOM displacement analysis, the ALD with retensioning and the ALD without retensioning were both found to have significantly higher displacement values than the FLD ($P < 0.002$). When comparing the retensioned ALD to the non-retensioned ALD, there was no significant difference ($P = 0.130$) [Figure 2].

In the ABM displacement analysis, there was a significant difference between the non-retensioned ALD and the FLD ($P = 0.018$). There was no significant difference between the retensioned ALD and the FLD ($P = 0.995$). When comparing the non-retensioned ALD to the retensioned ALD, there was no significant difference ($P = 0.317$) [Figure 3].

The secondary across-study meta-analysis of total displacement for the three groups (ALD-retensioned, ALD-not retensioned, and FLD) yielded similar results to those from the standardized mean differences. Those data are available online in the SI [Figure SI 1].

Meta-analysis failure load

The data extracted for the failure load analysis can be found on-line in the SI Data Source Tables ([Table SI 3a] for DOM failure load and [Table SI 3b] for ABM failure load).

In the DOM failure load analysis, the ALD with retensioning

Table 1: Lists the specific loading protocols for each study.

Study	Device only Model				Animal Model			
	Precycling	Loading Cycles	Loading Forces (low-high) (N)	Failure Testing	Precycling	Loading Cycles	Loading Forces (low-high) (N)	Failure Testing
Barrow et al	10 cycles 10-50N 1Hz	4500	10-250	20mm/min	No animal model			
Chang, M et al	10-75 N for 0.1Hz for 10 cycles	4500	100-400	50mm/min	10-75 N for 0.1Hz for 10 cycles	4500	100-400	50mm/min
Cheng, J et al	10 cycles 10-50N 1Hz	1000	50-250	50mm/hr	No animal model			
Eguchi et al	50N for 30 seconds	2000	50-250	60mm/min	50N for 30 sec	2000	50-250	60mm/min
Johnson et al	10-75 N for 0.1Hz for 10 cycles	1000	100-400	50mm/min	No animal model			
Noonan et al	10 cycles of 10N to -50N, then 25 cycles of 50-250N	4500	50-250	50mm/min	50N for 3sec, then 25 cycles of 50-250N	1000	50-250	50mm/min
Nye et al	No DOM study				25N tension for 30sec, then one pull from 50-250N	1000	50-250	50mm/min
Smith et al	No DOM study				10 cycles from 1mm to -2mm at 0.5Hz	1. 1000 cycles at 1mm to -2mm. 2. 1000 cycles at 10-250N		50mm/min
Petre et al	10 cycles 10-50N 0.1Hz	1000	50-250	50mm/min	10 cycles 10-50N 0.1Hz	1000	50-250	50mm/min
Gotschi et al	10 cycles 10-50N	5000	50-250	20mm/min	10 cycles 10-50N	1000	50-250	20mm/min
Singh et al	10cycles 20-70N 1Hz	5000	20-520	20mm/min	No animal model			
Conner et al	N/A				20 & 50 cycles 50-450N	2000	50-450	30mm/s
Glasbrenner et al	N/A				10 cycles 10-50N	2500	0,50,100,150,200,250	N/A
Rylander et al	N/A				N/A	1000	50-250	20mm/min
Ahmad et al	100N at 0.1mm/s	2000	50-500	0.2mm/s	N/A			
Kamelger et al	10 cycles 50-250	1000	50-250	20mm/min	10 cycles 50-250	1000	50-250	20mm/min

N/A means that it was not detailed in the study. mm: Millimeters, min: Minute, s: Second, Hz: Hertz, N: Newtons

and the ALD without retensioning both had significantly lower load to failure compared to the FLD ($P < 0.003$). There was no significant difference between the ALD without retensioning and the ALD with retensioning ($P = 0.716$) [Figure 4].

In the ABM failure load analysis, there was no significant difference between the ALD with retensioning and the FLD ($P = 0.500$). The ALD without retensioning had a significantly lower load to failure than the FLD ($P < 0.001$) and the ALD with retensioning ($P = 0.008$) [Figure 5].

The secondary across-study meta-analysis of failure forces for the three groups (ALD-retensioned, ALD-not retensioned, and FLD) yielded similar results to those from the standardized mean differences. Those data are available online in the SI [Figure SI 2].

DISCUSSION

This study demonstrated that isolated testing of both the ALD with retensioning and the ALD without retensioning had significantly higher cyclical displacement and significantly lower load to failure compared to the FLD. However, in the animal model data, there was no significant difference between the ALD with retensioning and the FLD. Meta-analyses regarding non-retensioned ALD demonstrated

the higher displacement and the lower load to failure in comparison to FLD. The results of this study suggest that the use of an ALD with intraoperative retensioning may produce a biomechanically equivalent construct to that of the FLD.

The outcomes of this analysis compare favorably and expand on the results of a recent meta-analysis.^[15] Houck *et al.* performed a meta-analysis comparing adjustable loops to fixed loops and found the adjustable loop had significantly higher displacement.^[15] In contrast to this analysis, Houck *et al.* found the ENDOBUTTON™ FLD to have a similar load to failure as the ToggleLoc™ adjustable loop but a significantly higher load to failure compared to the TightRope® RT adjustable loop.^[15] There have been concerns about the Houck study regarding the grouping of animal and DOMs as well as the exclusions of some biomechanical studies.^[27] In addition, their study did not consider the retensioning of the adjustable loop. The present study separates the comparison of animal and device only testing and includes a larger number of studies. This analysis demonstrates the high variability in the data even with this separation, with heterogeneity scores above 70% for all analyses. This may be due to differences in the devices themselves or the testing protocols (e.g., peak load level during the cyclic loading regime).




Figure 2: Forest plot for the meta-analysis results of standard mean displacement difference for the device only model studies. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “20” or “40” next to device refers to the length of the device loop in millimeters. RIGIDLOOP® A refers to adjustable RIGIDLOOP® device. RIGIDLOOP® NA refers to non-adjustable RIGIDLOOP® device. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group. CI: Confidence Interval




Figure 3: Forest plot for the meta-analysis of standardized mean displacement difference from animal model studies. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “15,” “20,” “21,” or “40” next to device refers to the length of the device loop in millimeters. “L” refers to placement on lateral cortex whereas “A” refers to placement on the anterior cortex. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group. CI: Confidence Interval.

Retensioning is an important property of the ALD and is recommended by multiple manufacturers.^[28,29] After the precycling period, the sutures can be pulled tight again to restore some of the displacement caused during precycling. An initial displacement of 0.75 mm after one pull at 250 N of force was reported by Petre *et al.*^[3] This displacement was 68% of the final lengthening seen after all the loading cycles were completed. When retensioning was done after the initial precycling, this reduced the initial length and, thus, the overall final displacement was significantly reduced.^[17] Noonan *et al.*,^[17] Nye *et al.*,^[18] and Smith *et al.*^[8] found that there was no significant difference in displacement between the FLD and the ALD that were retensioned. Conversely, Johnson *et al.*^[13] found that the ALD with retensioning displaced significantly more than the FLD. This study provides further evidence to support retensioning. In the current study, there was no significant difference between the FLD and the ALD with retensioning in the animal model analysis. In contrast, there was a significant difference between the FLD and ALD with retensioning in the DOM. The animal model is more likely to be clinically relevant as the porcine bone used in the included studies has been shown to closely replicate the biomechanical properties of the human knee.^[3,30,31] This analysis is also limited by the high heterogeneity of the studies. Tying the free suture ends of the adjustable loop is another strategy that may be used to limit the cyclical displacement.

There are two main adjustable-loop designs. One design uses a finger trap mechanism to prevent slippage once tension is applied while the other design uses a suture loop that locks once tension is applied. Once the button is retensioned after cyclical displacement, a knot can be tied to prevent the suture from sliding through the finger trap mechanism, or from slipping through the locking loop. This would create the equivalent of a fixed-loop construct. The evidence regarding the use of knotting is favorable but limited. Barrow *et al.*^[2] and Noonan *et al.*^[17] found that after securing the loose suture ends with knots that the adjustable loop displaced significantly less. Cheng *et al.*^[10] found no significant difference in displacement after knotting the ALD. In regard to failure load, Barrow *et al.*^[2] and Cheng *et al.*^[10] found that knotting made no difference, while Noonan *et al.*^[17] found that knotting significantly increased the failure load. This meta-analysis did not further investigate the utility of knotting the suture ends given the limited number of studies that studied this technique specifically. Despite the findings of this biomechanical analysis, there are data suggesting that ALD and FLD have similar clinical outcomes. ACL clinical failure is defined as absolute displacement >10 mm or >3 mm of side to side difference, which is reflected by total displacement measures.^[32] The cortical suspension device is supposed to maintain sufficient fixation so that graft tension is maintained until it is




Figure 4: Forest plot for standardized mean difference in failure load results from device only model studies. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “20” or “40” next to device refers to length of device loop in millimeters. RIGIDLOOP[®] A refers to adjustable RIGIDLOOP[®] device. RIGIDLOOP[®] NA refers to non-adjustable RIGIDLOOP[®] device. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group. CI: Confidence interval.




Figure 5: Forest plot of the meta-analysis results for standardized mean difference in failure load in animal model studies. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “15,” “20,” “21,” or “40” next to device refers to the length of the device loop in millimeters. “L” refers to placement on lateral cortex whereas “A” refers to placement on the anterior cortex. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group. CI: Confidence interval.

incorporated into native bone.^[2,33,34] If the suspension device lengthens during this initial period, then there would be a loss of graft tension which may also lead to tunnel widening and impaired healing.^[2,35] It should be cautioned that the additive effects of using dual suspensory fixation could lead to clinical failure. Smith *et al.*, the only biomechanical study to evaluate dual suspensory ALD, found that it displaced more than using a single ALD, but did not reach statistical significance.^[8] Supplementary tibial fixation, particularly in dual suspensory constructs, could be another method to potentially decrease cyclical elongation.^[36] Both the ALD and FLD demonstrated higher load to failure than what is thought to occur during rehabilitation. The forces that load the device during early rehabilitation and ambulation have been shown to be between 300 N and 590 N,^[6,14,33] while the grafts themselves have all been shown to be stronger than the native ACL with ultimate failure loads over 2300 N.^[37] Therefore, because all three groups showed failure loads above 590 N, it is probable that they are strong enough to withstand the stress of early rehabilitation. In support of this, two recent systematic reviews of clinical studies found no clinical differences between the retensioned ALD and the FLD.^[38,39] This discrepancy between biomechanical and clinical results could be partially explained by retensioning. Onggo *et al.*^[38] noted that all the ALD used in clinical studies were retensioned compared to few retensioned ALD in biomechanical studies. Other possible factors for the discrepancy between biomechanical and clinical studies include differing loading angles and the effect of the graft healing in the setting of early rehabilitation.^[38] Limitations of the present study include the experimental setup of studies. These animal studies typically used forces that were in line with the bone tunnels rather than the force vectors that might be acting on the knee *in vivo*. In biomechanical studies, the loading force is placed solely on the device. Conversely, in the native knee, the forces applied to the devices are not purely linear, with the graft-bone interface sharing the loading force.^[3,8,14] This means that the device most likely experiences lower forces *in vivo* and may not require as much strength as *in vitro*. The models lacked synovial fluid and tissue healing and were performed at time zero, which negates the effects of bone graft integration. The animal model consisted of fewer studies and could be underpowered to discover a significant difference. The heterogeneity of studies included is also a significant limitation. The I² values of our studies ranged from 71% to 90%. This meta-analysis combined data from studies comparing a range of ALD and FLD designs, assuming that these devices respond similarly which may not be the case. The studies also differed in their lengths of devices, cycling protocols, bone tunnel lengths, and graft lengths as well. These results may not translate clinically and require further clinical study.

CONCLUSION

Retensioning ALD improved their biomechanical performance in animal models but did not significantly improve their performance during device-only testing. However, there is significant heterogeneity in the pooled biomechanical results. Establishing a consistent testing protocol for these devices may help show whether retensioning of ALD alters the biomechanical response for these devices.

Acknowledgments

The authors gratefully acknowledge Theresa Atkinson, PhD for her assistance in the data analysis and manuscript preparation.

Declaration of patient consent

Patient's consent not required as there are no patients in this study.

Financial support and sponsorship

Nil.

Conflicts of interest

One or more of the authors has declared the following potential conflict of interest or source of funding: N.H.A. has received consulting fees from Pacira Pharmaceuticals, Smith and Nephew, Baudax Bio, AcelRx Pharmaceuticals, Trice Medical, Medical Device Business Services, Bioventus, KCI USA, and DePuy Synthes Products; royalties from Smith and Nephew; and hospitality payments from PAVmed and Acumed. J.N.L. has received support for education from Empire Medical, Smith and Nephew, and Arthrex and hospitality payments from Stryker Corporation and Smith and Nephew.

Use of artificial intelligence (AI)-assisted technology for manuscript preparation

The author(s) confirms that there was no use of artificial intelligence (AI)-assisted technology for assisting in the writing or editing of the manuscript and no images were manipulated using AI.

REFERENCES

1. Mall NA, Chalmers PN, Moric M, Tanaka MJ, Cole BJ, Bach BR Jr., *et al.* Incidence and trends of anterior cruciate ligament reconstruction in the United States. *Am J Sports Med* 2014;42:2363-70.
2. Barrow AE, Pilia M, Guda T, Kadrmas WR, Burns TC. Femoral suspension devices for anterior cruciate ligament reconstruction: Do adjustable loops lengthen? *Am J Sports Med* 2014;42:343-9.
3. Petre BM, Smith SD, Jansson KS, de Meijer PP, Hackett TR, LaPrade RF, *et al.* Femoral cortical suspension devices for soft tissue anterior cruciate ligament reconstruction: A comparative biomechanical study. *Am J Sports Med* 2013;41:416-22.
4. Walsh MP, Wijdicks CA, Parker JB, Hapa O, LaPrade RF. A comparison between a retrograde interference screw, suture button, and combined fixation on the tibial side in an all-inside anterior cruciate ligament reconstruction: A biomechanical study in a porcine model. *Am J Sports Med* 2009;37:160-7.
5. Stadelmaier DM, Lowe WR, Ilahi OA, Noble PC, Kohl HW 3rd. Cyclic pull-

- out strength of hamstring tendon graft fixation with soft tissue interference screws. Influence of screw length. *Am J Sports Med* 1999;27:778-83.
6. Kokkinakis M, Ashmore A, El-Guindi M. Intraoperative complications using the Bio-transfix femoral fixation implant in anterior cruciate ligament reconstruction. *Arch Orthop Trauma Surg* 2010;130:375-9.
 7. Ozer H, Ozgur A. Complications following hamstring anterior cruciate ligament reconstruction with femoral cross-pin fixation. *Arthroscopy* 2005;21:1407-8.
 8. Smith PA, Piepenbrink M, Smith SK, Bachmaier S, Bedi A, Wijdicks CA. Adjustable-versus fixed-loop devices for femoral fixation in ACL reconstruction: An *in vitro* full-construct biomechanical study of surgical technique-based tibial fixation and graft preparation. *Orthop J Sports Med* 2018;6:2325967118768743.
 9. Eguchi A, Ochi M, Adachi N, Deie M, Nakamae A, Usman MA. Mechanical properties of suspensory fixation devices for anterior cruciate ligament reconstruction: Comparison of the fixed-length loop device versus the adjustable-length loop device. *Knee* 2014;21:743-8.
 10. Cheng J, Paluvadi SV, Lee S, Yoo S, Song EK, Seon JK. Biomechanical comparisons of current suspensory fixation devices for anterior cruciate ligament reconstruction. *Int Orthop* 2018;42:1291-6.
 11. Conner CS, Perez BA, Morris RP, Buckner JW, Buford WL Jr, Ivey FM. Three femoral fixation devices for anterior cruciate ligament reconstruction: Comparison of fixation on the lateral cortex versus the anterior cortex. *Arthroscopy* 2010;26:796-807.
 12. Glasbrenner J, Domnick C, Raschke MJ, Willinghöfer T, Kittl C, Michel P, *et al.* Adjustable buttons for ACL graft cortical fixation partially fail with cyclic loading and unloading. *Knee Surg Sports Traumatol Arthrosc* 2019;27:2530-6.
 13. Johnson JS, Smith SD, LaPrade CM, Turnbull TL, LaPrade RF, Wijdicks CA. A biomechanical comparison of femoral cortical suspension devices for soft tissue anterior cruciate ligament reconstruction under high loads. *Am J Sports Med* 2015;43:154-60.
 14. Kamelger FS, Onder U, Schmoelz W, Tecklenburg K, Arora R, Fink C. Suspensory fixation of grafts in anterior cruciate ligament reconstruction: A biomechanical comparison of 3 implants. *Arthroscopy* 2009;25:767-76.
 15. Houck DA, Kraeutler MJ, McCarty EC, Bravman JT. Fixed-versus adjustable-loop femoral cortical suspension devices for anterior cruciate ligament reconstruction: A systematic review and meta-analysis of biomechanical studies. *Orthop J Sports Med* 2018;6:2325967118801762.
 16. Chang MJ, Bae TS, Moon YW, Ahn JH, Wang JH. A comparative biomechanical study of femoral cortical suspension devices for soft-tissue anterior cruciate ligament reconstruction: Adjustable-length loop versus fixed-length loop. *Arthroscopy* 2018;34:566-72.
 17. Noonan BC, Dines JS, Allen AA, Altchek DW, Bedi A. Biomechanical evaluation of an adjustable loop suspensory anterior cruciate ligament reconstruction fixation device: The value of retensioning and knot tying. *Arthroscopy* 2016;32:2050-9.
 18. Nye DD, Mitchell WR, Liu W, Ostrander RV. Biomechanical comparison of fixed-loop and adjustable-loop cortical suspensory devices for metaphyseal femoral-sided soft tissue graft fixation in anatomic anterior cruciate ligament reconstruction using a porcine model. *Arthroscopy* 2017;33:1225-32.e1.
 19. Ahmad SS, Hirschmann MT, Voumard B, Kohl S, Zysset P, Mukabeta T, *et al.* Adjustable loop ACL suspension devices demonstrate less reliability in terms of reproducibility and irreversible displacement. *Knee Surg Sports Traumatol Arthrosc* 2018;26:1392-8.
 20. Hozo SP, Djulbegovic B, Hozo I. Estimating the mean and variance from the median, range, and the size of a sample. *BMC Med Res Methodol* 2005;5:13.
 21. Rylander L, Brunelli J, Taylor M, Baldini T, Ellis B, Hawkins M, *et al.* A biomechanical comparison of anterior cruciate ligament suspensory fixation devices in a porcine cadaver model. *Clin Biomech (Bristol, Avon)* 2014;29:230-4.
 22. Viswanathan M, Ansari MT, Berkman ND, Chang S, Hartling L, McPheeters M, *et al.* Assessing the risk of bias of individual studies in systematic reviews of health care interventions. In: *Methods Guide for Effectiveness and Comparative Effectiveness Reviews*. Rockville, MD: Agency for Healthcare Research and Quality (US), AHRQ Methods for Effective Health Care; 2008.
 23. Higgins JP, Thompson SG. Quantifying heterogeneity in a meta-analysis. *Stat Med* 2002;21:1539-58.
 24. Higgins JP, Thompson SG, Deeks JJ, Altman DG. Measuring inconsistency in meta-analyses. *BMJ* 2003;327:557-60.
 25. Götschi T, Rosenberg G, Li X, Zhang C, Bachmann E, Snedeker JG, *et al.* Biomechanical evaluation of a novel loop retention mechanism for cortical graft fixation in ACL reconstruction. *Orthop J Sports Med* 2020;8:2325967120904322.
 26. Singh S, Ramos-Pascual S, Czerbak K, Malik M, Schranz PJ, Miles AW, *et al.* Biomechanical testing of fixed and adjustable femoral cortical suspension devices for ACL reconstruction under high loads and extended cyclic loading. *J Exp Orthop* 2020;7:27.
 27. Shang X, Li J, Li Q. Systematic review of fixed-versus adjustable-loop femoral cortical suspension devices for anterior cruciate ligament reconstruction: Letter to the editor. *Orthop J Sports Med* 2019;7:2325967119843504.
 28. GraftLink® ACL reconstruction with the ACL tightrope® II fixation: Surgical technique.
 29. ToggleLoc™ and toggleloc XL for ACL reconstruction using a quadruple semitendinosus graft.
 30. Seil R, Rupp S, Krauss PW, Benz A, Kohn DM. Comparison of initial fixation strength between biodegradable and metallic interference screws and a press-fit fixation technique in a porcine model. *Am J Sports Med* 1998;26:815-9.
 31. Xerogeanes JW, Fox RJ, Takeda Y, Kim HS, Ishibashi Y, Carlin GJ, *et al.* A functional comparison of animal anterior cruciate ligament models to the human anterior cruciate ligament. *Ann Biomed Eng* 1998;26:345-52.
 32. Alford JW, Bach BR Jr. Arthrometric aspects of anterior cruciate ligament surgery before and after reconstruction with patellar tendon grafts. *Tech Orthop* 2005;20:421-38.
 33. Kousa P, Järvinen TL, Vihavainen M, Kannus P, Järvinen M. The fixation strength of six hamstring tendon graft fixation devices in anterior cruciate ligament reconstruction. Part I: Femoral site. *Am J Sports Med* 2003;31:174-81.
 34. Shelton WR, Fagan BC. Autografts commonly used in anterior cruciate ligament reconstruction. *J Am Acad Orthop Surg* 2011;19:259-64.
 35. Rodeo SA, Kawamura S, Kim HJ, Dynybil C, Ying L. Tendon healing in a bone tunnel differs at the tunnel entrance versus the tunnel exit: An effect of graft-tunnel motion? *Am J Sports Med* 2006;34:1790-800.
 36. Eisen SH, Davidson PA, Rivenburgh DW. Supplemental tibial fixation for anterior cruciate ligament reconstruction. *Arthroscopy* 2008;24:1078-80.
 37. Dargel J, Gotter M, Mader K, Pennig D, Koebke J, Schmidt-Wiethoff R. Biomechanics of the anterior cruciate ligament and implications for surgical reconstruction. *Strategies Trauma Limb Reconstr* 2007;2:1-12.
 38. Onggo JR, Nambiar M, Pai V. Fixed-versus adjustable-loop devices for femoral fixation in anterior cruciate ligament reconstruction: A systematic review. *Arthroscopy* 2019;35:2484-98.
 39. Singh S, Shaunak S, Shaw SC, Anderson JL, Mandalia V. Adjustable loop femoral cortical suspension devices for anterior cruciate ligament reconstruction: A systematic review. *Indian J Orthop* 2020;54:426-43.

How to cite this article: Campbell ST, Liu J, Gowd AK, Amin N, Sardelli M, Morrison M. Does retensioning of adjustable-loop cortical suspension devices improve performance: A systematic review and meta-analysis. *J Arthrosc Surg Sports Med*. doi: 10.25259/JASSM_22_2022

SUPPLEMENTARY TABLES

Table SI 1: Rating for each question: minimal risk of bias 0; low risk of bias 1; moderate to high risk 2.

	<i>Is there a standardized sampling process for selecting the included devices?</i>	<i>Is testing consistent with all tests performed the same way?</i>	<i>Is the sample size calculation described justified?</i>	<i>Is there a well described, consistent method used to analyze the test data to produce the outcome measures?</i>	<i>Are all test results reported? If an outlier data point is removed is this well described and justified?</i>	RISK OF BIAS #	RATING
<i>Ahmed</i>	1	0	1	0	0	2	low
<i>Chang M</i>	0	0	1	0	0	1	low
<i>Cheng J</i>	1	0	1	0	0	2	low
<i>Connor</i>	0	1	1	0	0	2	low
<i>Eguchi</i>	2	0	2	0	0	2	moderate
<i>Glassbrenner</i>	0	0	1	0	0	1	low
<i>Gotschi</i>	1	0	0	0	0	1	low
<i>Johnson</i>	2	0	1	0	0	3	low
<i>Kamelger</i>	1	0	1	0	0	2	low
<i>Noonan</i>	1	0	1	0	2	4	moderate
<i>Nye</i>	1	0	1	0	1	3	moderate
<i>Petre</i>	2	0	1	0	0	3	low
<i>Rylander</i>	0	0	1	0	0	1	low
<i>Singh</i>	1	0	1	0	0	2	low
<i>Smith</i>	1	1	1	0	0	3	low

Table SI 2a: Raw data for displacement in DOM.

ALD	DOM displacement								
	N	Failure load (N)	SD	Retensioned or not	FLD/Control	N	Failure load (N)	SD	Retensioned or not
Johnson Tightrope NRT	8	2.2	0.62	NRT	Johnson Rigidloop	8	1.09	0.16	N/A
Johnson ToggleLoc NRT	8	3.69	2.39	NRT	Johnson XO	8	1.65	0.43	N/A
Johnson Tightrope RT	8	1.81	0.51	RT	Johnson Endobutton	8	1.05	0.05	N/A
Johnson ToggleLoc RT	8	3.22	1.41	RT	Johnson Endobutton	8	1.05	0.05	N/A
Noonan Tightrope NRT	5	0.96	0.07	NRT	Noonan Endobutton	5	0.42	0.08	N/A
Noonan Tightrope NRT+K	5	0.63	0.14	NRT	Noonan Endobutton	5	0.42	0.08	N/A
Noonan Tightrope RT	5	0.73	0.10	RT	Noonan Endobutton	5	0.42	0.08	N/A
Noonan Tightrope RT+K	5	0.38	0.09	RT	Noonan Endobutton	5	0.42	0.08	N/A
Noonan Tightrope NRT Unloaded	5	4.22	2.68	NRT	Noonan Endobutton	5	0.42	0.08	N/A
Noonan Tightrope RT Unloaded	5	0.51	0.11	RT	Noonan Endobutton	5	0.42	0.08	N/A
Cheng J Tightrope NRT	8	1.56	0.08	NRT	Cheng, J Endobutton	8	0.76	0.06	N/A
Cheng J Tightrope NRT+K	8	1.38	0.24	NRT	Cheng, J Endobutton	8	0.76	0.06	N/A
Cheng, J Graft Max NRT	8	2.11	0.57	NRT	Cheng, J Endobutton	8	0.76	0.06	N/A
Chang, M Tightrope NRT	6	1.99	0.4	NRT	Chang, M Endobutton	6	0.79	0.05	N/A
Eguchi Tightrope NRT	10	4.05	1.16	NRT	Eguchi Endobutton	10	2.03	0.31	N/A
Petre Tightrope NRT	5	1.1	0.2	NRT	Petre XO	5	1.2	0.17	N/A
Petre ToggleLoc NRT	5	2.18	0.31	NRT	Petre Endobutton	5	0.42	0.08	N/A
Gotschi Tightrope NRT	11	1.26	0.52	NRT	Gotschi FlippTack	8	1.823	0.17	N/A
Gotschi Variloop NRT	11	0.85	0.08	NRT	Gotschi FlippTack	8	1.823	0.17	N/A
Singh Ultrabutton RT	5	2.66	0.28	RT	Singh G-Lok	5	1.46	0.25	N/A
Singh Rigidloop RT	5	1.51	0.16	RT	Singh G-Lok	5	1.46	0.25	N/A
Singh ProCinch RT	5	1.6	0.09	RT	Singh G-Lok	5	1.46	0.25	N/A
Kamelger ToggleLoc 20 NRT	6	0.66	0.12	NRT	Kamelger Endobutton 20	6	0.15	0.01	N/A
Kamelger ToggleLoc 40 NRT	6	0.76	0.06	NRT	Kamelger Endobutton 40	6	0.22	0.03	N/A
Ahmed Rigidloop A NRT	6	1.67	0.27	NRT	Ahmed Endobutton	6	1.07	0.06	N/A
Ahmed ProCinch NRT	6	3.57	2.05	NRT	Ahmed Retrobutton	6	0.69	0.04	N/A
Ahmed Ultrabutton NRT	6	3.14	0.66	NRT	Ahmed Rigidloop NA	6	1.22	0.09	N/A

NRT: Non-retensioned, RT: Retensioned, N/A: Not applicable, SD: Standard deviation, DOM: Device-only model, ALD: Adjustable-loop device, FLD: Fixed-loop device, K: Knotted. N: Number of devices. Number "20" or "40" next to device refers to the length of the device loop in millimeters.

Table SI 2b: Raw data for displacement in ABM.

ALD	Animal model displacement								
	N	Failure load (N)	SD	Retensioned or not	FLD/Control	N	Failure load (N)	SD	Retensioned or not
Nye Tightrope RT	10	5.09	0.87	RT	Nye Endobutton	10	5.07	0.56	N/A
Nye Toggleloc RT	10	7.44	1.63	RT	Nye Endobutton	10	5.07	0.56	N/A
Smith Tightrope RT	8	2.78	0.85	RT	Smith Retrobutton	8	2.85	1.03	N/A
Smith Ultrabutton RT	8	2.76	0.45	RT	Smith Endobutton	8	2.85	0.74	N/A
Smith Graft Max RT	8	4.13	1.46	RT	Smith Endobutton	8	2.85	0.74	N/A
Noonan Tightrope NRT	5	2.7	0.5	NRT	Noonan Endobutton	5	3.00	0.30	N/A
Noonan Tightrope RT+K	5	1.5	0.30	RT	Noonan Endobutton	5	3.00	0.30	N/A
Chang, M Tightrope NRT	6	15.65	2.43	NRT	Chang, M Endobutton	6	14.88	1.79	N/A
Eguchi Tightrope 21 NRT	10	7.74	2.52	NRT	Eguchi Endobutton	10	5.88	1.06	N/A
Eguchi Tightrope 15 NRT	10	6.39	2.32	NRT	Eguchi Endobutton	10	5.88	1.06	N/A
Petre Tightrope NRT	10	4.47	0.65	NRT	Petre XO Button	10	3.5	0.5	N/A
Petre Toggleloc NRT	10	6.02	1.9	NRT	Petre Endobutton	10	3.37	0.27	N/A
Glasbrenner Ultrabutton NRT	8	8.1	1.50	NRT	Glasbrenner Endobutton	8	4.4	0.30	N/A
Glasbrenner Tightrope NRT	8	6.1	1.40	NRT	Glasbrenner FlippTack	8	4.1	0.60	N/A
Glasbrenner Graft Max NRT	8	4.7	1.00	NRT	Glasbrenner Endobutton	8	4.4	0.30	N/A
Conner ToggleLoc L NRT	2	6.52	0.72	NRT	Conner Endobutton L	3	4.84	0.72	N/A
Conner ToggleLoc A NRT	5	5.46	0.95	NRT	Conner Endobutton A	5	3.55	0.57	N/A
Gotschi Tighrope NRT	8	3.921	0.68	NRT	Gotschi FlippTack	8	4.666	0.72	N/A
Gotschi Variloop NRT	8	2.183	0.32	NRT	Gotschi FlippTack	8	4.666	0.72	N/A
Kamelger ToggleLoc 20 NRT	8	1.56	0.44	NRT	Kamelger Endobutton 20	8	0.66	0.19	N/A
Kamelger ToggleLoc 40 NRT	8	1.37	0.33	NRT	Kamelger Endobutton 40	8	0.5	0.13	N/A

NRT: Non-retensioned, RT: Retensioned, N/A: Not applicable, SD: Standard deviation, ALD: Adjustable-loop device, FLD: Fixed loop device, ABM: Animal bone model, K: Knotted. N: Number of devices. Number "15," "20," "21," or "40" next to device refers to the length of the device loop in millimeters.

Table SI 3a: Raw data for failure load in DOM.

ALD	DOM failure load								
	N	Failure load (N)	SD	Retensioned or not	FLD/Control	N	Failure load (N)	SD	Retensioned or not
Johnson Tightrope NRT	8	784	45	NRT	Johnson Rigidloop	8	1976	229	N/A
Johnson ToggleLoc NRT	8	1995	217	NRT	Johnson XO	8	2218	114	N/A
Johnson Tightrope RT	8	1020	421	RT	Johnson Endobutton	8	1530	180	N/A
Johnson ToggleLoc RT	8	2231	511	RT	Johnson Endobutton	8	1530	180	N/A
Noonan Tightrope NRT	5	886	39	NRT	Noonan Endobutton	5	1384.0	71.0	N/A
Noonan Tightrope RT	5	815	35	RT	Noonan Endobutton	5	1384.0	71.0	N/A
Noonan Tightrope NRT+K	5	1038	71	NRT	Noonan Endobutton	5	1384.0	71.0	N/A
Noonan Tightrope RT+K	5	1205	35	RT	Noonan Endobutton	5	1384.0	71.0	N/A
Noonan Tightrope NRT Unloaded	5	833	151	NRT	Noonan Endobutton	5	1384.0	71.0	N/A
Noonan Tightrope RT Unloaded	5	1057	156	RT	Noonan Endobutton	5	1384.0	71.0	N/A
Cheng J Tightrope NRT	8	800.9	112.5	NRT	Cheng, J Endobutton	8	1204.0	127.0	N/A
Cheng J Tightrope NRT+K	8	868	51.9	NRT	Cheng, J Endobutton	8	1204.0	127.0	N/A
Cheng, J Graft Max NRT	8	914.2	36.9	NRT	Cheng, J Endobutton	8	1204.0	127.0	N/A
Chang, M Tightrope NRT	6	925	38.12	NRT	Chang, M Endobutton	6	1410.0	118.4	N/A
Eguchi Tightrope NRT	10	866	53	NRT	Eguchi Endobutton	10	1430.0	148.0	N/A
Petre Tightrope NRT	5	841	55	NRT	Petre XO	5	2230	252	N/A
Petre ToggleLoc NRT	5	1561	112	NRT	Petre Endobutton	5	1456	130	N/A
Gotschi Tightrope NRT	11	827.31	34.02	NRT	Gotschi FlippTack	8	1317.12	120.99	N/A
Gotschi Variloop NRT	11	817.78	96.52	NRT	Gotschi FlippTack	8	1317.12	120.99	N/A
Singh Ultrabutton RT	5	1903	81.00	RT	Singh G-Lok	5	2178	118.00	N/A
Singh Rigidloop RT	5	1835	179.00	RT	Singh G-Lok	5	2178	118.00	N/A
Singh ProCinch RT	5	1456	137.00	RT	Singh G-Lok	5	2178	118.00	N/A
Kamelger ToggleLoc 20 NRT	6	1557.6	97.50	NRT	Kamelger Endobutton 20	6	1074.6	119.40	N/A
Kamelger ToggleLoc 40 NRT	6	1454.2	170.00	NRT	Kamelger Endobutton 40	6	1202.2	59.7	N/A
Ahmed Rigidloop A NRT	6	1676.62	58.13	NRT	Ahmed Endobutton	6	1524.11	135.13	N/A
Ahmed ProCinch NRT	6	811.32	25.42	NRT	Ahmed Retrobutton	6	1259.2	89.74	N/A
Ahmed Ultrabutton NRT	6	1710.69	58.10	NRT	Ahmed Rigidloop NA	6	1931.34	140.93	N/A

NRT: Non-retensioned, RT: Retensioned, N/A: Not applicable, SD: Standard deviation, DOM: Device-only model, ALD: Adjustable-loop device, FLD: Fixed-loop device, K: Knotted. N: Number of devices. Number "20" or "40" next to device refers to the length of the device loop in millimeters.

Table SI 3b: Raw data for failure load in ABM.

ALD	N	Failure Load (N)	SD	Animal model failure load					
				Retensioned or not	FLD/Control	N	Failure Load (N)	SD	Retensioned or not
Nye Tightrope RT	10	801.1	56.3	RT	Nye Endobutton	10	803.9	92.2	N/A
Nye ToggleLoc RT	10	682.1	182.4	RT	Nye Endobutton	10	803.9	92.2	N/A
Smith Tightrope RT	8	958	40	RT	Smith Retrobutton	8	689	134	N/A
Smith Ultrabutton RT	8	746	180	RT	Smith Endobutton	8	712.0	78.0	N/A
Smith Graft Max RT	8	761	150	RT	Smith Endobutton	8	712.0	78.0	N/A
Noonan Tightrope NRT	5	786	166	NRT	Noonan Endobutton	5	866.0	102.0	N/A
Noonan Tightrope RT+K	5	818	168	RT	Noonan Endobutton	5	866.0	102.0	N/A
Chang, M Tightrope NRT	6	888	90.31	NRT	Chang, M Endobutton	6	843.0	111.4	N/A
Eguchi Tightrope 21 NRT	10	880	60	NRT	Eguchi Endobutton	10	1115.0	274.0	N/A
Eguchi Tightrope 15 NRT	10	860	70	NRT	Eguchi Endobutton	10	1115.0	274.0	N/A
Petre Tightrope NRT	10	859	43	NRT	Petre XO Button	10	1748	140	N/A
Petre ToggleLoc NRT	10	1334	81	NRT	Petre Endobutton	10	1456	101	N/A
Conner ToggleLoc L NRT	2	876	207.00	NRT	Conner Endobutton L	3	987	305	N/A
Conner ToggleLoc A NRT	5	913	82.00	NRT	Conner Endobutton A	5	1191	150.00	N/A
Gotschi Tightrope NRT	8	873.47	56.34	NRT	Gotschi Flipptack	8	1201.98	234.04	N/A
Gotschi Variloop NRT	8	849.49	68.25	NRT	Gotschi Flipptack	8	1201.98	234.04	N/A
Kamelger ToggleLoc 20 NRT	8	943.4	199.00	NRT	Kamelger Endobutton 20	8	1024.7	75.00	N/A
Kamelger ToggleLoc 40 NRT	8	945.4	233.10	NRT	Kamelger Endobutton 40	8	1122.6	95.3	N/A
Rylander ToggleLoc NRT	9	559.7	101.30	NRT	Rylander Endobutton	10	716.7	128.20	N/A

NRT: Non-retensioned, RT: Retensioned, N/A: Not applicable, SD: Standard deviation, ALD: Adjustable-loop device, ABM: Animal bone models. Number "15," "20," "21," or "40" next to device refers to the length of the device loop in millimeters.

SUPPLEMENTARY FIGURES

Secondary meta-analysis


Animal combined displacement

Figure SI 1: Forest plot for secondary analysis of total displacement in animal models. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “15,” “20,” “21,” or “40” next to the device refers to the length of the device loop in millimeters. “L” refers to placement on the lateral cortex, whereas “A” refers to placement on the anterior cortex. CI: Confidence interval. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group.

Combined animal failure load




Figure SI 2: Forest plot for secondary analysis of load to failure in animal models. NRT: Not retensioned, RT: Retensioned, K: Knotted. Unloaded refers to the Noonan *et al.* protocol that featured smaller lower limit forces during cyclical testing.^[17] Number “15,” “20,” “21,” or “40” next to the device refers to the length of the device loop in millimeters. “L” refers to placement on the lateral cortex, whereas “A” refers to placement on the anterior cortex. CI: Confidence interval. Yellow is the subgroup standardized mean difference. Blue is combined standardized mean difference. Black squares refer to the mean for that study and group.